

Product datasheet for **GA110752**

Human AP1AR activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	AP1AR CRISPRa kit - CRISPR gene activation of human adaptor related protein complex 1 associated regulatory protein
Symbol:	AP1AR
Components:	GA110752G1, AP1AR gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110752G2, AP1AR gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110752G3, AP1AR gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q63HQ0</u>
Synonyms:	2C18; C4orf16; gamma-BAR; GBAR; PRO0971
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)