


Product datasheet for **GA110754**

Human CAMK2N1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CAMK2N1 CRISPRa kit - CRISPR gene activation of human calcium/calmodulin dependent protein kinase II inhibitor 1
Symbol:	CAMK2N1
Components:	GA110754G1, CAMK2N1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110754G2, CAMK2N1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110754G3, CAMK2N1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7Z7J9</u>
Synonyms:	PRO1489
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)