

Product datasheet for **GA110769**

Human DDX43 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DDX43 CRISPRa kit - CRISPR gene activation of human DEAD-box helicase 43
Symbol:	DDX43
Components:	GA110769G1, DDX43 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110769G2, DDX43 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110769G3, DDX43 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NXZ2</u>
Synonyms:	CT13; HAGE
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)