

Product datasheet for **GA110789**

Human Plexin A3 (PLXNA3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLXNA3 CRISPRa kit - CRISPR gene activation of human plexin A3
Symbol:	PLXNA3
UniProt ID:	P51805
Synonyms:	6.3; HSSEXGENE; PLXN3; PLXN4; XAP-6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)