

Product datasheet for **GA110793**

Human DNAH3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DNAH3 CRISPRa kit - CRISPR gene activation of human dynein axonemal heavy chain 3
Symbol:	DNAH3
UniProt ID:	Q8TD57
Synonyms:	DLP3; dnahc3-b; DNAHC3B; HEL-36; Hsadhc3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)