

Product datasheet for **GA110800**

Human **FAM48A (SUPT20H)** activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SUPT20H CRISPRa kit - CRISPR gene activation of human SPT20 homolog, SAGA complex component
Symbol:	SUPT20H
UniProt ID:	<u>Q8NEM7</u>
Synonyms:	C13; C13orf19; FAM48A; FP757; P38IP; SPT20
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)