


Product datasheet for **GA110813**

Human CDKN2AIP activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CDKN2AIP CRISPRa kit - CRISPR gene activation of human CDKN2A interacting protein
Symbol:	CDKN2AIP
Components:	GA110813G1, CDKN2AIP gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110813G2, CDKN2AIP gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110813G3, CDKN2AIP gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NXV6</u>
Synonyms:	CARF
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)