

Product datasheet for **GA110814**

Human **TENT5A** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TENT5A CRISPRa kit - CRISPR gene activation of human terminal nucleotidyltransferase 5A
Symbol:	TENT5A
Components:	GA110814G1, TENT5A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110814G2, TENT5A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110814G3, TENT5A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96IP4</u>
Synonyms:	C6orf37; XTP11
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)