

## Product datasheet for **GA110826**

### Human ASAP3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	ASAP3 CRISPRa kit - CRISPR gene activation of human ArfGAP with SH3 domain, ankyrin repeat and PH domain 3
Symbol:	ASAP3
UniProt ID:	<a href="#">Q8TDY4</a>
Synonyms:	ACAP4; CENTB6; DDEFL1; UPLC1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)