

Product datasheet for **GA110837**

Human ZNF407 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF407 CRISPRa kit - CRISPR gene activation of human zinc finger protein 407
Symbol:	ZNF407
Components:	GA110837G1, ZNF407 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110837G2, ZNF407 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110837G3, ZNF407 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9C0G0
Synonyms:	FLJ13839; FLJ20307; KIAA1703
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)