

Product datasheet for **GA110840**

Human **LRRC40** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	LRRC40 CRISPRa kit - CRISPR gene activation of human leucine rich repeat containing 40
Symbol:	LRRC40
Components:	GA110840G1, LRRC40 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110840G2, LRRC40 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110840G3, LRRC40 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H9A6</u>
Synonyms:	dj677H15.1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

