

Product datasheet for **GA110841**

Human KIAA1333 (G2E3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	G2E3 CRISPRa kit - CRISPR gene activation of human G2/M-phase specific E3 ubiquitin protein ligase
Symbol:	G2E3
UniProt ID:	Q7L622
Synonyms:	KIAA1333; PHF7B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)