

## Product datasheet for **GA110873**

### Human Smek1 (PPP4R3A) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PPP4R3A CRISPRa kit - CRISPR gene activation of human protein phosphatase 4 regulatory subunit 3A
Symbol:	PPP4R3A
UniProt ID:	<u><a href="#">Q6IN85</a></u>
Synonyms:	FLFL1; KIAA2010; MSTP033; PP4R3; PP4R3A; SMEK1; smk-1; smk1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)