

Product datasheet for **GA110882**

Human DSU (MREG) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MREG CRISPRa kit - CRISPR gene activation of human melanoregulin
Symbol:	MREG
Components:	GA110882G1, DSU gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110882G2, DSU gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110882G3, DSU gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N565</u>
Synonyms:	DSU; WDT2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)