


Product datasheet for **GA110911**

Human TSR1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSR1 CRISPRa kit - CRISPR gene activation of human TSR1 ribosome maturation factor
Symbol:	TSR1
Components:	GA110911G1, TSR1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110911G2, TSR1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110911G3, TSR1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q2NL82</u>
Synonyms:	FLJ10534; KIAA1401; MGC131829
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)