

Product datasheet for **GA110912**

Human IQCC activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IQCC CRISPRa kit - CRISPR gene activation of human IQ motif containing C
Symbol:	IQCC
Components:	GA110912G1, IQCC gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110912G2, IQCC gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110912G3, IQCC gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q4KMZ1</u>
Synonyms:	RP4-622L5.6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)