

Product datasheet for **GA110940**

Human Transmembrane protein 30A (TMEM30A) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | TMEM30A CRISPRa kit - CRISPR gene activation of human transmembrane protein 30A |
| Symbol: | TMEM30A |
| UniProt ID: | Q9NV96 |
| Synonyms: | C6orf67; CDC50A |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)