

## Product datasheet for **GA110962**

### Human WDR52 (CFAP44) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	CFAP44 CRISPRa kit - CRISPR gene activation of human cilia and flagella associated protein 44
Symbol:	CFAP44
Components:	GA110962G1, WDR52 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110962G2, WDR52 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110962G3, WDR52 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u><a href="#">Q96MT7</a></u>
Synonyms:	WDR52
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)