

Product datasheet for **GA110963**

Human C6orf70 (ERMARD) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | ERMARD CRISPRa kit - CRISPR gene activation of human ER membrane associated RNA degradation |
| Symbol: | ERMARD |
| UniProt ID: | <u>Q5T6L9</u> |
| Synonyms: | C6orf70; dj266L20.3; PVNH6 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)