

Product datasheet for **GA110973**

Human C1orf103 (LRIF1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LRIF1 CRISPRa kit - CRISPR gene activation of human ligand dependent nuclear receptor interacting factor 1
Symbol:	LRIF1
UniProt ID:	Q5T3J3
Synonyms:	C1orf103; RIF1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)