

Product datasheet for **GA110976**

Human PCID2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PCID2 CRISPRa kit - CRISPR gene activation of human PCI domain containing 2
Symbol:	PCID2
Components:	GA110976G1, PCID2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110976G2, PCID2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110976G3, PCID2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5JVF3</u>
Synonyms:	F10
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)