

## Product datasheet for **GA110979**

### Human CACNA2D3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	CACNA2D3 CRISPRa kit - CRISPR gene activation of human calcium voltage-gated channel auxiliary subunit alpha2delta 3
Symbol:	CACNA2D3
UniProt ID:	<a href="#">Q8IZS8</a>
Synonyms:	HSA272268
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)