

Product datasheet for **GA110980**

Human **SCN3B** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SCN3B CRISPRa kit - CRISPR gene activation of human sodium voltage-gated channel beta subunit 3
Symbol:	SCN3B
UniProt ID:	<u>Q9NY72</u>
Synonyms:	ATFB16; BRGDA7; HSA243396; SCN3B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)