

Product datasheet for **GA110993**

Human TSNAXIP1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSNAXIP1 CRISPRa kit - CRISPR gene activation of human translin associated factor X interacting protein 1
Symbol:	TSNAXIP1
Components:	GA110993G1, TSNAXIP1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110993G2, TSNAXIP1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110993G3, TSNAXIP1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q2TAA8</u>
Synonyms:	TXI1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)