

Product datasheet for **GA111001**

Human NRIP (DCAF6) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DCAF6 CRISPRa kit - CRISPR gene activation of human DDB1 and CUL4 associated factor 6
Symbol:	DCAF6
Components:	GA111001G1, NRIP gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111001G2, NRIP gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111001G3, NRIP gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q58WW2
Synonyms:	1200006M05Rik; ARCAP; IQWD1; MSTP055; NRIP; PC326
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)