

Product datasheet for **GA111004**

Human TMEM111 (EMC3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EMC3 CRISPRa kit - CRISPR gene activation of human ER membrane protein complex subunit 3
Symbol:	EMC3
UniProt ID:	Q9P0I2
Synonyms:	POB; TMEM111
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)