

Product datasheet for **GA111008**

Human EAPP activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EAPP CRISPRa kit - CRISPR gene activation of human E2F associated phosphoprotein
Symbol:	EAPP
Components:	GA111008G1, EAPP gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111008G2, EAPP gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111008G3, EAPP gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q56P03
Synonyms:	BM036; C14orf11
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)