

Product datasheet for **GA111009**

Human CENPN activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CENPN CRISPRa kit - CRISPR gene activation of human centromere protein N
Symbol:	CENPN
Components:	GA111009G1, CENPN gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111009G2, CENPN gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111009G3, CENPN gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96H22
Synonyms:	BM039; C16orf60; CENP-N; ICEN32
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)