

## Product datasheet for **GA111011**

### Human WWC3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	WWC3 CRISPRa kit - CRISPR gene activation of human WWC family member 3
Symbol:	WWC3
Components:	GA111011G1, WWC3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111011G2, WWC3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111011G3, WWC3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u><a href="#">Q9ULE0</a></u>
Synonyms:	BM042
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)