

Product datasheet for **GA111017**

Human C9ORF46 (PLGRKT) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLGRKT CRISPRa kit - CRISPR gene activation of human plasminogen receptor with a C-terminal lysine
Symbol:	PLGRKT
UniProt ID:	Q9HBL7
Synonyms:	AD025; C9orf46; MDS030; Plg-R(KT); PLG-RKT
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)