

Product datasheet for **GA111019**

Human PEN2 (PSENEN) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSENEN CRISPRa kit - CRISPR gene activation of human presenilin enhancer, gamma-secretase subunit
Symbol:	PSENEN
UniProt ID:	<u>Q9NZ42</u>
Synonyms:	MDS033; MSTP064; PEN-2; PEN2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)