

Product datasheet for **GA111020**

Human **TEX2** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TEX2 CRISPRa kit - CRISPR gene activation of human testis expressed 2
Symbol:	TEX2
Components:	GA111020G1, TEX2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111020G2, TEX2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111020G3, TEX2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8IWB9</u>
Synonyms:	HT008; TMEM96
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)