

Product datasheet for **GA111024**

Human TMEM165 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM165 CRISPRa kit - CRISPR gene activation of human transmembrane protein 165
Symbol:	TMEM165
Components:	GA111024G1, TMEM165 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111024G2, TMEM165 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111024G3, TMEM165 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9HC07</u>
Synonyms:	CDG2K; FT27; GDT1; TMPT27; TPARL
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)