

Product datasheet for **GA111034**

Human PBK activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PBK CRISPRa kit - CRISPR gene activation of human PDZ binding kinase
Symbol:	PBK
UniProt ID:	Q96KB5
Synonyms:	CT84; HEL164; Nori-3; SPK; TOPK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)