

Product datasheet for **GA111037**

Human **WSB2** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	WSB2 CRISPRa kit - CRISPR gene activation of human WD repeat and SOCS box containing 2
Symbol:	WSB2
Components:	GA111037G1, WSB2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111037G2, WSB2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111037G3, WSB2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NYS7</u>
Synonyms:	SBA2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)