

Product datasheet for **GA111039**

Human ZNF167 (ZKSCAN7) activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | ZKSCAN7 CRISPRa kit - CRISPR gene activation of human zinc finger with KRAB and SCAN domains 7 |
| Symbol: | ZKSCAN7 |
| UniProt ID: | <u>Q9P0L1</u> |
| Synonyms: | ZFP; ZNF64; ZNF167; ZNF448; ZSCAN39 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)