

Product datasheet for **GA111047**

Human **UNC45A** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	UNC45A CRISPRa kit - CRISPR gene activation of human unc-45 myosin chaperone A
Symbol:	UNC45A
UniProt ID:	Q9H3U1
Synonyms:	GC-UNC45; GCUNC-45; GCUNC45; IRO039700; SMAP-1; SMAP1; UNC-45A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)