


Product datasheet for **GA111049**

Human THSD1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	THSD1 CRISPRa kit - CRISPR gene activation of human thrombospondin type 1 domain containing 1
Symbol:	THSD1
Components:	GA111049G1, THSD1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111049G2, THSD1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111049G3, THSD1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NS62</u>
Synonyms:	TMTSP; UNQ3010
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)