

Product datasheet for **GA111069**

Human LIN37 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LIN37 CRISPRa kit - CRISPR gene activation of human lin-37 DREAM MuvB core complex component
Symbol:	LIN37
Components:	GA111069G1, LIN37 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111069G2, LIN37 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111069G3, LIN37 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96GY3</u>
Synonyms:	F25965; lin-37; ZK418.4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)