

Product datasheet for **GA111079**

Human **SLC25A40** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SLC25A40 CRISPRa kit – CRISPR gene activation of human solute carrier family 25 member 40
Symbol:	SLC25A40
Components:	GA111079G1, SLC25A40 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111079G2, SLC25A40 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111079G3, SLC25A40 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8TBP6</u>
Synonyms:	MCFP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

