

Product datasheet for **GA111081**

Human RAG1AP1 (SLC50A1) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | SLC50A1 CRISPRa kit - CRISPR gene activation of human solute carrier family 50 member 1 |
| Symbol: | SLC50A1 |
| UniProt ID: | Q9BRV3 |
| Synonyms: | HsSWEET1; RAG1AP1; SCP; slv; SWEET1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)