


Product datasheet for **GA111094**

Human TMEM234 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM234 CRISPRa kit - CRISPR gene activation of human transmembrane protein 234
Symbol:	TMEM234
Components:	GA111094G1, TMEM234 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111094G2, TMEM234 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111094G3, TMEM234 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8WY98</u>
Synonyms:	AASL548; C1orf91; dj622L5.7; PRO1105; RP4-622L5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)