

Product datasheet for **GA111096**

Human PCDHGC4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PCDHGC4 CRISPRa kit – CRISPR gene activation of human protocadherin gamma subfamily C, 4
Symbol:	PCDHGC4
Components:	GA111096G1, PCDHGC4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111096G2, PCDHGC4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111096G3, PCDHGC4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y5F7</u>
Synonyms:	PCDH-GAMMA-C4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

