

## Product datasheet for **GA111100**

### Human PCDHGB3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PCDHGB3 CRISPRa kit - CRISPR gene activation of human protocadherin gamma subfamily B, 3
Symbol:	PCDHGB3
UniProt ID:	<a href="#">Q9Y5G1</a>
Synonyms:	PCDH-GAMMA-B3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:

