


Product datasheet for **GA111109**

Human PCDHGA4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PCDHGA4 CRISPRa kit - CRISPR gene activation of human protocadherin gamma subfamily A, 4
Symbol:	PCDHGA4
Components:	GA111109G1, PCDHGA4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111109G2, PCDHGA4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111109G3, PCDHGA4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y5G9</u>
Synonyms:	PCDH-GAMMA-A4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)