

Product datasheet for **GA111113**

Human PCDHB15 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | PCDHB15 CRISPRa kit - CRISPR gene activation of human protocadherin beta 15 |
| Symbol: | PCDHB15 |
| UniProt ID: | Q9Y5E8 |
| Synonyms: | PCDH-BETA15 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)