

Product datasheet for **GA111114**

Human PCDHB14 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PCDHB14 CRISPRa kit - CRISPR gene activation of human protocadherin beta 14
Symbol:	PCDHB14
Components:	GA111114G1, PCDHB14 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111114G2, PCDHB14 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111114G3, PCDHB14 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9Y5E9
Synonyms:	PCDH-BETA14
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)