

Product datasheet for **GA111117**

Human PCDHB11 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PCDHB11 CRISPRa kit – CRISPR gene activation of human protocadherin beta 11
Symbol:	PCDHB11
Components:	GA111117G1, PCDHB11 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111117G2, PCDHB11 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111117G3, PCDHB11 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9Y5F2
Synonyms:	ME2; PCDH-BETA11
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

