

Product datasheet for **GA111119**

Human PCDHB9 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PCDHB9 CRISPRa kit - CRISPR gene activation of human protocadherin beta 9
Symbol:	PCDHB9
Components:	GA111119G1, PCDHB9 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111119G2, PCDHB9 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111119G3, PCDHB9 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y5E1</u>
Synonyms:	PCDH-BETA9; PCDH3H
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)