

Product datasheet for **GA111124**

Human Protocadherin beta 3 (PCDHB3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PCDHB3 CRISPRa kit - CRISPR gene activation of human protocadherin beta 3
Symbol:	PCDHB3
Components:	GA111124G1, Protocadherin beta 3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111124G2, Protocadherin beta 3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111124G3, Protocadherin beta 3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9Y5E6
Synonyms:	PCDH-BETA3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)